



11128/FA

On board medical suction unit with autoclavable canister



FEATURES

- Rescue vehicles On board medical suction unit
- Maximum aspiration flow 30lt/min.
- Maximum aspiration Power 800 mbar.
- 1000 ml Autoclavabile canister with disposable (model LN)
- 1000 ml autoclavable canister with safety valve and protection filter directly fixed onto the cap (model FA)
- 12 Volt vehicle Connection Cable
- Dimension cm. 9x8,5x6 depth.
- In compliance with all existing norms.

SEE ALSO



11128/LN

On board medical suction unit with
disposable liner
