



# 7080/4RG Proof Mercury Lite

---

Ambulance cot in aluminium with variable height and four swivel castors



---

## FEATURES

Our "Mercury Lite" stretcher is available in the 4 swivel castors version, ensuring excellent manoeuvrability. This device is ideal for narrow spaces. The two front castors can be unblocked by means of an intuitive mechanism. They can revert to their original non-rotating position. The fastening of the stretcher inside of an emergency vehicle is ensured by our certified product 955 Proof (see "[Fixing Systems](#)"). The stretcher is supplied with mattress, straps and object holder.



Device certified according to the European harmonized safety standards UNI EN 1865 – UNI EN 1789.



Device certified according to the European harmonized safety standards UNI EN 1865 - UNI EN 1789.

\*All Pictures shown are for illustration purpose only.

MAIN FEATURES

10G Certification	
Material	Painted aluminium HDPE (High-density polyethylene)
Multilevel	
Lying surface	With Trendelenburg
Intermediate height (cm)	72
Max height (cm)	95
Max loading height (cm)	66
Min loading height (cm)	64
Folded height (cm)	31
Width (cm)	57
Length (cm)	197
Length with folded trolley (cm)	165
Weight (kg)	40
Loading capacity (kg)	250
Type	Barella autocaricante
Wheel type	4 Big Rotating Wheels
Usage	Transport
Size	Adults

IMAGES



SEE ALSO



**Mercury 7070/4RG PROOF**  
Stainless steel variable height  
stretcher certified with 4 rotating  
wheels



**Mercury Lite 7082/4RG Proof**  
Certified aluminium variable height  
stretcher high version with 4  
rotating wheels



**Mercury 7072/4RG PROOF**  
Certified stainless steel variable  
height stretcher high version with 4  
rotating wheels



**Mercury Cinque 7075/4RG PROOF**  
Certified 5 levels self-loading  
stretcher with 4 swivel wheels



**Mercury Cinque 7077/4RG PROOF**  
Certified 5 levels self-loading  
stretcher with 4 swivel wheels -  
high version



**Mercury EL 7092 PROOF**  
Certified self-loading stretcher



**Mercury EL 7094 PROOF**  
Certified self-loading stretcher -  
high version

## ACCESSORIES



**12080**  
"MPB" wide bariatric platform for self-loading stretchers Frog Mercury and Gecko complete of belts with reel and fixing for stretcher



**12080-001**  
Monting kit for bariatric platform on self loading stretcher



**12128**  
Adjustable system for headrest of Meber-X and Mercury stretchers



**12154**  
"BEAR PODS" Bariatric containment system for self-loading stretchers



**12018N**  
Telescopic foldable phleboclysis stand for self loading stretchers



**955 PROOF**  
EN 1789 fixing system for Proof roll-in stretchers



**12070/BG**  
"COVER OX 2-3 BG"- professional cylinder holder for 2 and 3 lt oxygen bottles



**12070/BR**  
COVER OX2-3 BR bariatric cylinder holder with belt and orange MeBer logo for 2 and 3 lt oxygen bottles



**7013/A**  
3 pcs. Skay fire-proof electrowelded mattress



**7014/A**  
4 pcs. Skay fire-proof electrowelded mattress



**12100/N**  
Black electrowelded shaped cushion with MeBer logo for electrowelded mattresses



**12106**  
Mattress for AIR-COM stretchers - Anti-bedsore mattress with alternate segments



**939**  
Instrument support for stretchers



**938**  
Oximeb oxygen cylinder holder



**12054**  
Oximeb Maxi - Oxygen cylinder holder



**12084-001**  
Cover for bariatric wings set



**12170**  
Black mattress for right bariatric wing for self-loading stretchers



**12172**  
Black mattress for left bariatric wing for self-loading stretchers



**707**  
Red 6-handle fireproof cover with MeBer tape for self-loading stretcher



**608/A-MEB**  
2 pieces orange belt with metallic hook and MeBer band (L. 155 cm)



**608/N-MEB**  
2 pieces black belt with metallic hook and MeBer band (L. 155 cm)



**612/AN-2/MEB**



**612/N-2/MEB**



**12025**  
Black multi usage object holder for automatic loading stretcher

Adjustable thoracic abdominal belt  
orange/black with MeBer band and  
4 points fastening for self loading  
stretchers and sedan chairs

Adjustable thoracic abdominal belt  
black with MeBer band with 4  
points fastening for self loading  
stretchers and sedan chairs




**12072/BG**

"COVER OX 5 BG"- professional  
cylinder holder for 5 lt oxygen  
bottles



**12072/BR**

"COVER OX 5 BR"- professional  
cylinder holder for 5 lt oxygen  
bottles



**12200**

Bariatric wings set for Mercury self-  
loading stretchers complete with  
mattresses, belts with reel and  
assembly kit

## STANDARD ACCESSORIES



**7013/N**

3 pcs. Skay fire-proof electrowelded  
mattress



**608/N-MEB**

2 pieces black belt with metallic  
hook and MeBer band (L. 155 cm)



**612/N-2/MEB**

Adjustable thoracic abdominal belt  
black with MeBer band with 4  
points fastening for self loading  
stretchers and sedan chairs



**12025**

Black multi usage object holder for  
automatic loading stretcher