



7082/4RG Proof Mercury Lite

Certified aluminium variable height stretcher high version with 4 rotating wheels



FEATURES


4 rotating wheels version is available for "Mercury Lite" which grants extremely maneuverability. This product is ideal for narrow spaces. Operating an easy command the two frontal wheels can be unlocked or reset to their original configuration. Fixing of stretcher inside the rescue vehicle is granted by certified system Art. 955 Proof (see "[Fixing Systems](#)"). Supplied complete with mattress and belts. This version has the same features of model 7080/4RG Proof, but it allows to load on ambulances with higher loading height.



Device certified according to the European harmonized safety standards UNI EN 1865 - UNI EN 1789.

*All Pictures shown are for illustration purpose only.

MAIN FEATURES

10G Certification	
Material	Painted aluminium
Multilevel	
Lying surface	With Trendelenburg
Intermediate height (cm)	75
Max height (cm)	99
Max loading height (cm)	71
Min loading height (cm)	69
Folded height (cm)	31
Width (cm)	57
Length (cm)	197
Length with folded trolley (cm)	165
Weight (kg)	40,5
Loading capacity (kg)	250

IMAGES



SEE ALSO



Mercury 7070/4RG PROOF
Stainless steel variable height
stretcher certified with 4 rotating
wheels



Mercury Lite 7080/4RG Proof
Ambulance cot in aluminium with
variable height and four swivel
castors



Mercury EL 7092 PROOF
Certified self-loading stretcher



Mercury EL 7094 PROOF
Certified self-loading stretcher -
high version

ACCESSORIES



938
Oximeb oxygen cylinder holder



939
Instrument support for stretchers



955 PROOF
EN 1789 fixing system for Proof roll-in stretchers



7014/A
4 pcs. Skay fire-proof electrowelded mattress



12054
Oximeb Maxi - Oxygen cylinder holder



608/N-MEB
2 pieces black belt with metallic hook and MeBer band (L. 155 cm)



612/AN-2/MEB
Adjustable thoracic abdominal belt orange/black with MeBer band and 4 points fastening for self loading stretchers and sedan chairs



12080
"MPB" wide bariatric platform for self-loading stretchers Frog Mercury and Gecko complete of belts with reel and fixing for stretcher



12100/N
Black electrowelded shaped cushion with MeBer logo for electowelded mattresses



12106
Mattress for AIR-COM stretchers - Anti-bedsores mattress with alternate segments



12080-001
Monting kit for bariatric platform on self loading stretcher



12018N
Telescopic foldable phlebotomy stand for self loading stretchers



12154
"BEAR PODS" Bariatric containment system for self-loading stretchers



12128
Adjustable system for headrest of Meber-X and Mercury stretchers



12070/BR
COVER OX2-3 BR professional cylinder holder with belt and orange MeBer logo for 2 and 3 lt oxygen bottles



12070/BG
"COVER OX 2-3 BG"- professional cylinder holder for 2 and 3 lt oxygen bottles



12072/BR



12072/BG



12025
Black multi usage object holder for automatic loading stretcher



707
Red 6-handle fireproof cover with MeBer tape for self-loading stretcher

"COVER OX 5 BR"- professional cylinder holder for 5 lt oxygen bottles

"COVER OX 5 BC"- professional cylinder holder for 5 lt oxygen bottles



12170

Black mattress for right bariatric wing for self-loading stretchers



12172

Black mattress for left bariatric wing for self-loading stretchers



608/A-MEB

2 pieces orange belt with metallic hook and MeBer band (L. 155 cm)




612/N-2/MEB

Adjustable thoracic abdominal belt black with MeBer band with 4 points fastening for self loading stretchers and sedan chairs



12084-001

Cover for bariatric wings set



12200

Bariatric wings set for Mercury self-loading stretchers complete with mattresses, belts with reel and assembly kit

STANDARD ACCESSORIES



612/N-2/MEB

Adjustable thoracic abdominal belt black with MeBer band with 4 points fastening for self loading stretchers and sedan chairs



7013/N

3 pcs. Skay fire-proof electrowelded mattress



12025

Black multi usage object holder for automatic loading stretcher



608/N-MEB

2 pieces black belt with metallic hook and MeBer band (L. 155 cm)