



7082/4RG Proof Mercury Lite

Certified aluminium variable height stretcher high version with 4 rotating wheels



FEATURES



4 rotating wheels version is available for "Mercury Lite" which grants extremely maneuverability. This product is ideal for narrow spaces. Operating an easy command the two frontal wheels can be unlocked or reset to their original configuration. Fixing of stretcher inside the rescue vehicle is granted by certified system Art. 955 Proof (see "[Fixing Systems](#)"). Supplied complete with mattress and belts. This version has the same features of model 7080/4RG Proof, but it allows to load on ambulances with higher loading height.



Device certified according to the European harmonized safety standards UNI EN 1865 - UNI EN 1789.

*All Pictures shown are for illustration purpose only.

MAIN FEATURES

10G Certification	
Material	Painted aluminium HDPE (High-density polyethylene)
Multilevel	
Lying surface	With Trendelenburg
Intermediate height (cm)	75
Max height (cm)	99
Max loading height (cm)	71
Min loading height (cm)	69
Folded height (cm)	31
Width (cm)	57
Length (cm)	197
Length with folded trolley (cm)	165
Weight (kg)	40,5
Loading capacity (kg)	250

IMAGES



SEE ALSO



Mercury Lite 7080/4RG Proof
Ambulance cot in aluminium with
variable height and four swivel
castors



Mercury 7070/4RG PROOF
Stainless steel variable height
stretcher certified with 4 rotating
wheels



Mercury EL 7092 PROOF
Certified self-loading stretcher



Mercury EL 7094 PROOF
Certified self-loading stretcher -
high version

ACCESSORIES



12072/BG
"COVER OX 5 BG"- professional
cylinder holder for 5 lt oxygen
bottles



12106
Mattress for AIR-COM stretchers -
Anti-bedsore mattress with
alternate segments



707
Red 6-handle fireproof cover with
MeBer tape for self-loading
stretcher



938
Oximeb oxygen cylinder holder



12054
Oximeb Maxi - Oxygen cylinder
holder



12084-001
Cover for bariatric wings set



12170
Black mattress for right bariatric
wing for self-loading stretchers



12172
Black mattress for left bariatric
wing for self-loading stretchers



7013/A
3 pcs. Skay fire-proof electrowelded
mattress



7014/A
4 pcs. Skay fire-proof electrowelded
mattress



12100/N
Black electrowelded shaped
cushion with MeBer logo for
electrowelded mattresses



955 PROOF
EN 1789 fixing system for Proof roll-
in stretchers



608/A-MEB
2 pieces orange belt with metallic
hook and MeBer band (L. 155 cm)



608/N-MEB
2 pieces black belt with metallic
hook and MeBer band (L. 155 cm)



612/AN-2/MEB
Adjustable thoracic abdominal belt
orange/black with MeBer band and
4 points fastening for self loading
stretchers and sedan chairs



612/N-2/MEB
Adjustable thoracic abdominal belt
black with MeBer band with 4
points fastening for self loading
stretchers and sedan chairs



12025
Black multi usage object holder for
automatic loading stretcher



12080
"MPB" wide bariatric platform for
self-loading stretchers Frog
Mercury and Gecko complete of
belts with reel and fixing for
stretcher



12080-001
Monting kit for bariatric platform
on self loading stretcher



12154
"BEAR PODS" Bariatric
containment system for self-
loading stretchers



12018N

Telescopic foldable phlebotomy stand for self-loading stretchers



939

Instrument support for stretchers



12128

Adjustable system for headrest of Meber-X and Mercury stretchers



12070/BG

"COVER OX 2-3 BG"- professional cylinder holder for 2 and 3 lt oxygen bottles




12070/BR

COVER OX2-3 BR professional cylinder holder with belt and orange MeBer logo for 2 and 3 lt oxygen bottles



12072/BR

"COVER OX 5 BR"- professional cylinder holder for 5 lt oxygen bottles



12200

Bariatric wings set for Mercury self-loading stretchers complete with mattresses, belts with reel and assembly kit

STANDARD ACCESSORIES



12025

Black multi usage object holder for automatic loading stretcher



7013/N

3 pcs. Skay fire-proof electrowelded mattress



612/N-2/MEB

Adjustable thoracic abdominal belt black with MeBer band with 4 points fastening for self-loading stretchers and sedan chairs



608/N-MEB

2 pieces black belt with metallic hook and MeBer band (L. 155 cm)