



9216 Gemini

Cervical collar set 6 pieces with carrying case



FEATURES

GEMINI cervical collar is a highly professional instrument indispensable for the emergency operator. Composed by 2 parts, it is positioned in the front and rear part of the neck allowing thus a perfect immobilization of the rachis. Central hole for tracheotomy, easy to be used and adjust thanks to the rip-off fastening. Made with stamped polyethylene coloured to identify various size and E.V.A. washable inner lining fixed by plastic rivets. Available in 6 different size.

MAIN FEATURES

Trauma



IMAGES



SEE ALSO



Gemini 9200
Rigid cervical collar in 2 pieces size
INFANT



Gemini 9202
Rigid cervical collar in 2 pieces size
PEDIATRIC



Gemini 9204
Rigid cervical collar in 2 pieces size
NECKLESS



Gemini 9206
Rigid cervical collar in 2 pieces size
SHORT



Gemini 9208
Rigid cervical collar in 2 pieces size
REGULAR



Gemini 9210
Rigid cervical collar in 2 pieces size
TALL
