



## 9216 Gemini

---

Cervical collar set 6 pieces with carrying case



---

### FEATURES

GEMINI cervical collar is an highly professional instrument indispensable for the emergency operator. Composed by 2 parts, it is positioned in the front and rear part of the neck allowing thus a perfect immobilization of the rachis. Central hole for tracheotomy, easy to be used and adjust thanks to the rip-off fastening. Made with stamped polyethylene coloured to identify various size and E.V.A. washable inner lining fixed by plastic rivets. Available in 6 different size.

---

### MAIN FEATURES

---

Trauma



IMAGES



SEE ALSO



**Gemini 9200**  
Rigid cervical collar in 2 pieces size  
INFANT



**Gemini 9202**  
Rigid cervical collar in 2 pieces size  
PEDIATRIC



**Gemini 9204**  
Rigid cervical collar in 2 pieces size  
NECKLESS



**Gemini 9206**  
Rigid cervical collar in 2 pieces size  
SHORT



**Gemini 9208**  
Rigid cervical collar in 2 pieces size  
REGULAR



**Gemini 9210**  
Rigid cervical collar in 2 pieces size  
TALL