



## 9200 Gemini

---

Rigid cervical collar in 2 pieces size INFANT



---

### FEATURES

GEMINI cervical collar is an highly professional instrument indispensable for the emergency operator. Composed by 2 parts, it is positioned in the front and rear part of the neck allowing thus a perfect immobilization of the rachis. Central hole for tracheotomy, easy to be used and adjust thanks to the rip-off fastening. Made with stamped polyethylene coloured to identify various size and E.V.A. washable inner lining fixed by plastic rivets. Available in 6 different size.

---

### MAIN FEATURES

Trauma	<input checked="" type="checkbox"/>
Pediatric	<input checked="" type="checkbox"/>
Max length (cm)	24,5
Max height (cm)	17,5

---

### SEE ALSO



**Gemini 9216**  
Cervical collar set 6 pieces with  
carrying case

---