



9202 Gemini

Rigid cervical collar in 2 pieces size PEDIATRIC



FEATURES

GEMINI cervical collar is an highly professional instrument indispensable for the emergency operator. Composed by 2 parts, it is positioned in the front and rear part of the neck allowing thus a perfect immobilization of the rachis. Central hole for tracheotomy, easy to be used and adjust thanks to the rip-off fastening. Made with stamped polyethylene coloured to identify various size and E.V.A. washable inner lining fixed by plastic rivets. Available in 6 different size.

MAIN FEATURES

Trauma	<input checked="" type="checkbox"/>
Pediatric	<input checked="" type="checkbox"/>
Max length (cm)	25
Length (cm)	24,5

SEE ALSO



Gemini 9216
Cervical collar set 6 pieces with carrying case
