



9206 Gemini


Rigid cervical collar in 2 pieces size SHORT



FEATURES

GEMINI cervical collar is a highly professional instrument indispensable for the emergency operator. Composed by 2 parts, it is positioned in the front and rear part of the neck allowing thus a perfect immobilization of the rachis. Central hole for tracheotomy, easy to be used and adjust thanks to the rip-off fastening. Made with stamped polyethylene coloured to identify various size and E.V.A. washable inner lining fixed by plastic rivets. Available in 6 different size.

MAIN FEATURES

Trauma	
Max length (cm)	29
Length (cm)	28
Height (cm)	15
Max height (cm)	19

SEE ALSO



Gemini 9213

Cervical collar set 3 pieces with carrying case Art. 9204 - 9206 - 9208



Gemini 9214

Cervical collar set 4 pieces with carrying case Art. 9204 - 9206 - 9208 - 9210



Gemini 9216

Cervical collar set 6 pieces with carrying case