



## 9208 Gemini


Rigid cervical collar in 2 pieces size REGULAR



### FEATURES

GEMINI cervical collar is an highly professional instrument indispensable for the emergency operator. Composed by 2 parts, it is positioned in the front and rear part of the neck allowing thus a perfect immobilization of the rachis. Central hole for tracheotomy, easy to be used and adjust thanks to the rip-off fastening. Made with stamped polyethylene coloured to identify various size and E.V.A. washable inner lining fixed by plastic rivets. Available in 6 different size.

### MAIN FEATURES

|                 |   |
|-----------------|---|
| Trauma          |  |
| Length (cm)     | 31  |
| Height (cm)     | 16  |
| Max height (cm) | 20  |

### SEE ALSO



#### Gemini 9213

Cervical collar set 3 pieces with carrying case Art. 9204 - 9206 - 9208



#### Gemini 9214

Cervical collar set 4 pieces with carrying case Art. 9204 - 9206 - 9208 - 9210



#### Gemini 9216

Cervical collar set 6 pieces with carrying case