



## 9210 Gemini


Rigid cervical collar in 2 pieces size TALL



### FEATURES

GEMINI cervical collar is an highly professional instrument indispensable for the emergency operator. Composed by 2 parts, it is positioned in the front and rear part of the neck allowing thus a perfect immobilization of the rachis. Central hole for tracheotomy, easy to be used and adjust thanks to the rip-off fastening. Made with stamped polyethylene coloured to identify various size and E.V.A. washable inner lining fixed by plastic rivets. Available in 6 different size.

### MAIN FEATURES

Trauma	
Max length (cm)	36
Length (cm)	35
Height (cm)	19

### SEE ALSO



#### Gemini 9214

Cervical collar set 4 pieces with carrying case Art. 9204 - 9206 - 9208 - 9210



#### Gemini 9216

Cervical collar set 6 pieces with carrying case