



9236 Artis

Cervical collar set 6 pieces with carrying case



FEATURES

Rigid cervical collar in 1 piece made with stamped polyethylene coloured to identify various sizes and E.V.A. washable inner lining fixed by plastic rivets. Adjustable rip-off fastening and central hole for tracheotomy. Available in 6 different sizes.

MAIN FEATURES

Trauma



SEE ALSO



Artis 9220
Rigid cervical collar in 1 piece size
INFANT



Artis 9222
Rigid cervical collar in 1 piece size
PEDIATRIC



Artis 9224
Rigid cervical collar in 1 piece size
NECKLESS



Artis 9226
Rigid cervical collar in 1 piece size
SHORT



Artis 9228
Rigid cervical collar in 1 piece size
REGULAR



Artis 9230
Rigid cervical collar in 1 piece size
TALL