



9233 Artis

Cervical collar set 3 pieces with carrying case
Art. 9224 - 9226 - 9228



FEATURES

Rigid cervical collar in 1 piece made with stamped polyethylene coloured to identify various sizes and E.V.A. washable inner lining fixed by plastic rivets. Adjustable rip-off fastening and central hole for tracheotomy. Available in 6 different sizes.
Art. 9224 - 9226 - 9228

MAIN FEATURES

Trauma



SEE ALSO



Artis 9224
Rigid cervical collar in 1 piece size
NECKLESS



Artis 9226
Rigid cervical collar in 1 piece size
SHORT



Artis 9228
Rigid cervical collar in 1 piece size
REGULAR