



9222 Artis

Rigid cervical collar in 1 piece size PEDIATRIC



FEATURES

Rigid cervical collar in 1 piece made with stamped polyethylene coloured to identify various sizes and E.V.A. washable inner lining fixed by plastic rivets. Adjustable rip-off fastening and central hole for tracheotomy. Available in 6 different sizes.

MAIN FEATURES

Trauma	
Pediatric	

SEE ALSO



Artis 9234

Cervical collar set 4 pieces with carrying case. Art. 9220 - 9224 - 9226 - 9228



Artis 9236

Cervical collar set 6 pieces with carrying case