



## 9224 Artis

Rigid cervical collar in 1 piece size NECKLESS



### FEATURES

Rigid cervical collar in 1 piece made with stamped polyethylene coloured to identify various sizes and E.V.A. washable inner lining fixed by plastic rivets. Adjustable rip-off fastening and central hole for tracheotomy. Available in 6 different sizes.

### MAIN FEATURES

Trauma



### SEE ALSO



#### Artis 9233

Cervical collar set 3 pieces with carrying case Art. 9224 - 9226 - 9228



#### Artis 9234

Cervical collar set 4 pieces with carrying case. Art. 9220 - 9224 - 9226 - 9228



#### Artis 9234

Cervical collar set 4 pieces with carrying case. Art. 9220 - 9224 - 9226 - 9228



#### Artis 9236

Cervical collar set 6 pieces with carrying case