



9130 - ME.BER. BABY FIX


Fixing system with belts for transport of children inside the ambulance



FEATURES

Me.Ber. Baby Fix is a device for transport of children, adaptable to all model of stretcher by mean of 6 belts. Suitable for newborn and children up to 130 cm and 35 kg, allows a safe transport thanks to the restraint system at 5 points. Supplied with cover art. 12034. **Certified UNI EN 1789.** Length 74 x Head width 39 x Feet width 20 cm

MAIN FEATURES

Pediatric	
Length (cm)	74
Width (cm)	39
Weight (kg)	1,3
Usage	Transport

SEE ALSO



9128

Me.Ber. Vacuum Baby Fix

Vacuum mattress for fixing system
with belts for transport of children
inside the ambulance

ACCESSORIES



12034

Red empty carrying bag for "Baby
Fix" and "Vacuum Baby Fix"
