



630/N Maxima

Black painted anti trauma scoop stretcher



FEATURES

Keeping the very same features of Art.630, thanks to a special paint (designed for optimal durability), this stretcher can be used to differentiate operators or a specific area of operation.



Device certified according to the European harmonized safety standards UNI EN 1865 - UNI EN 1789.

- Effective hooking and unhooking device to open/close the stretcher.
- Adjustable length with automatic locking device.
- Bag with n. 3 quick unhooking safety belts to immobilize the patient (603/MEB).
- Easy arrangement inside the ambulance when folded.
- Completely built with light-weight and sturdy aluminium.

MAIN FEATURES

Size	Adults
Usage	Rescue
Height (cm)	7
Width (cm)	42
Length (cm)	167
Max length (cm)	202
Folded length (cm)	120
Epaisseur pliée (cm)	9
Weight (kg)	10
Loading capacity (kg)	170
10G Certification	
Material	Painted aluminium

SEE ALSO



Maxima 630
Grey anodized anti trauma scoop
stretcher



Maxima 630/A
Orange painted anti trauma scoop
stretcher



Maxima 630/B
Blue painted anti trauma scoop
stretcher



Maxima 630/G
Yellow painted anti trauma scoop
stretcher



Maxima 630/N
Black painted anti trauma scoop
stretcher



Maxima 630/N
Black painted anti trauma scoop
stretcher



Maxima 630/V
Green painted anti trauma scoop
stretcher



Maxima 631/G
Yellow painted anti trauma scoop
stretcher with head/trunk in one
piece



Maxima 631
Grey anodized anti trauma scoop
stretcher with head/trunk in 1 pc

ACCESSORIES



616
Certified vertical fastening system
for scoop stretchers



616/A
Certified vertical fixing system for
scoop stretchers



614
Certified horizontal fastening
system for scoop stretchers



630/CS
Carrying case for scoop stretchers



Double 624
Head immobilizer for scoop
stretcher



603/MEB
Nylon bag containing 3 belts Art.
607/Ameb