



## 631/G Maxima

Yellow painted anti trauma scoop  
stretcher with head/trunk in one piece



### FEATURES

It keeps same features of Art. 631. This stretcher helps differentiating operators of a specific area of operation thanks to its special yellow paint (designed for optimal durability).



Device certified according to the European harmonized safety standards UNI EN 1865 - UNI EN 1789

### MAIN FEATURES

Size	Adults
Usage	Rescue
Height (cm)	7
Width (cm)	42
Length (cm)	167
Max length (cm)	202
Folded length (cm)	120
Epaisseur pliée (cm)	9
Weight (kg)	9,5
Loading capacity (kg)	170
Type	
10G Certification	
Material	Painted aluminium

SEE ALSO



**Maxima 631**  
Grey anodized anti trauma scoop  
stretcher with head/trunk in 1 pc

ACCESSORIES



**Double 624**  
Head immobilizer for scoop  
stretcher



**614**  
Certified horizontal fastening  
system for scoop stretchers



**616**  
Certified vertical fastening system  
for scoop stretchers



**616/A**  
Certified vertical fixing system for  
scoop stretchers



**630/CS**  
Carrying case for scoop stretchers



**614**  
Certified horizontal fastening  
system for scoop stretchers



**603/MEB**  
Nylon bag containing 3 belts Art.  
607/Ameb