



7013/N

3 pcs. Skay fire-proof electrowelded mattress



FEATURES

Mattress for self loading stretcher made with Skay and fireproof padding. Totally electrowelded and water proof, can be used on almost all MeBer self loading stretchers with Trendelenburg including Cronos (not suitable for stretcher Evox and stretcher 1515).

Dimensions: cm 186 x 47 x 9.

MAIN FEATURES

Length (cm)	186
Width (cm)	47
Thickness (cm)	9
Material	Skay

SEE ALSO



mattress



^{*}All Pictures shown are for illustration purpose only.

ACCESSORY OF



MBS-01 MEBER BARIATRIC STRETCHER 7270 Self loading stretcher for bariatric transport



Mercury Cinque 7075/4RG PROOF Certified 5 levels self-loading stretcher with 4 swivel wheels



Mercury Cinque 7077/4RG PROOF Certified 5 levels self-loading stretcher with 4 swivel wheels high version



Mercury Lite Cinque 7095/4RG Proof Certified aluminium 5 levels selfloading stretcher with 4 swivel wheels



Mercury Lite Cinque 7095/4RG Proof Certified aluminium 5 levels selfloading stretcher with 4 swivel wheels



Mercury Lite Cinque 7097/4RG Proof Certified aluminium 5 levels selfloading stretcher with 4 swivel wheels - high version



Mercury EL 7092 PROOF Certified self-loading stretcher



Mercury EL 7092 PROOF Certified self-loading stretcher



Mercury EL 7094 PROOF Certified self-loading stretcher high version



7265 Proof Aluminium self loading stretcher "Frog lite Advance" with exctractable handles and Trendelenburg



Mercury Lite 7082/4RG Proof Certified aluminium variable height stretcher high version with 4 rotating wheels