



7013/N

3 pcs. Skay fire-proof electrowelded mattress



FEATURES

Mattress for self loading stretcher made with Skay and fireproof padding. Totally electrowelded and water proof, **can be used on almost all MeBer self loading stretchers with Trendelenburg including Cronos (not suitable for stretcher Evox and stretcher 1515).**

Dimensions: cm 186 x 47 x 9.

*All Pictures shown are for illustration purpose only.

MAIN FEATURES

Length (cm)	186
Width (cm)	47
Thickness (cm)	9
Material	Skay

SEE ALSO



7013/A
3 pcs. Skay fire-proof electrowelded mattress



7013/A
3 pcs. Skay fire-proof electrowelded mattress

ACCESSORY OF



MBS-01 MEBER BARIATRIC STRETCHER 7270
Self loading stretcher for bariatric transport



Mercury Cinque 7075/4RG PROOF
Certified 5 levels self-loading stretcher with 4 swivel wheels



Mercury Cinque 7077/4RG PROOF
Certified 5 levels self-loading stretcher with 4 swivel wheels - high version



Mercury Lite Cinque 7095/4RG Proof
Certified aluminium 5 levels self-loading stretcher with 4 swivel wheels



Mercury Lite Cinque 7095/4RG Proof
Certified aluminium 5 levels self-loading stretcher with 4 swivel wheels



Mercury Lite Cinque 7097/4RG Proof
Certified aluminium 5 levels self-loading stretcher with 4 swivel wheels - high version



Mercury EL 7092 PROOF
Certified self-loading stretcher



Mercury EL 7092 PROOF
Certified self-loading stretcher



Mercury EL 7094 PROOF
Certified self-loading stretcher - high version



7265 Proof
Aluminium self loading stretcher "Frog lite Advance" with extractable handles and Trendelenburg



Mercury Lite 7082/4RG Proof
Certified aluminium variable height stretcher high version with 4 rotating wheels