



## 910/GIN Winner

Self loading trolley for incubators



### FEATURES

Designed to load and transport incubators on rescue vehicles. This trolley is part of "Winner" range, keeping the all same features. Made of anodized aluminium.

### MAIN FEATURES

Pediatric	
Height (cm)	87,5
Max loading height (cm)	65
Min loading height (cm)	45
Folded height (cm)	26
Width (cm)	56
Length (cm)	194
Length with folded trolley (cm)	162
Weight (kg)	28
Material	Anodized aluminium
Loading capacity (kg)	170

---

## SEE ALSO

---



**Winner 910/GIN PROOF**  
Self loading trolley for incubators  
predisposed for certified fixing  
system

---

---

## ACCESSORIES

---



**Rear fixing system**  
Rear fixing system for roll-in H=148  
mm



**Front blocking device**  
Front blocking device for automatic  
loading stretcher



**12048**  
Certified interface for incubator  
fixation for self loading trolleys



**12050**  
Extended trolley incubator  
interface

---