




9040 Galaxy-P


Yellow spinal board with pins



FEATURES

Same features as “Galaxy” spinal boards Art. 9010, 9012 and 9015 equipped with 10 plastic PINS useful for a fast fixing of restraint systems with karabiners. Supplied with 3 belts (605/MEB).

MAIN FEATURES

Length (cm)	183
Width (cm)	41
Thickness (cm)	5
Loading capacity (kg)	150
Weight (kg)	5,8
Trauma	
Usage	Immobilization
Material	ABS

ACCESSORIES



Jelly 629

Orange standard head immobilizer
for spinal board



574

Blue bag with 3 belts art. 572

STANDARD ACCESSORIES



605/MEB

50 mm orange belt in one piece
with plastic buckle and ME.BER.
ribbon