



9128 - Me.Ber. Vacuum Baby Fix



Vacuum mattress for fixing system with belts for transport of children inside the ambulance



FEATURES

Me.Ber. Vacuum Baby Fix is a vacuum mattress to be used with Art. 9130. Extreme lightness and sturdiness, security belts and perfect encircling of patient make Art. 9128 unique. **Certified UNI EN 1789**, X-ray transparent, nonmagnetic.
Length 79 x Width 45 cm

MAIN FEATURES

Trauma	
Pediatric	
Length (cm)	79
Width (cm)	45
Weight (kg)	0,75
Usage	Transport

IMAGES



SEE ALSO



9130

ME.BER. BABY FIX

Fixing system with belts for
transport of children inside the
ambulance

ACCESSORIES



12034

Red empty carrying bag for "Baby
Fix" and "Vacuum Baby Fix"
