



853 Antares

Cervical collar set 3 pieces with carrying case
Art. 844 - 846 - 848



FEATURES

ANTARES cervical collar is an highly professional instrument indispensable for the emergency operator. It is composed by 2 pieces to be put on the front and back part of the neck enabling thus the perfect immobilization of the rachis. Central hole for tracheotomy and easy positioning thanks to the rip-off fastening which different colours indicate the various sizes. It is made with transparent polyethylene with E.V.A. washable inner lining fixed by plastic rivets. Available in 6 different sizes.
Art. 844 - 846 - 848


MAIN FEATURES

Usage	Immobilization
Trauma	

SEE ALSO



Sirius 874
Cervical collar set 4 pieces with
carrying case Art. 862 - 864 - 866 -
868



Sirius 876
Cervical collar set 6 pieces with
carrying case



Antares 844
Rigid cervical collar in 2 pieces size
STOUT



Antares 846
Rigid cervical collar in 2 pieces size
STANDARD



Antares 848
Rigid cervical collar in 2 pieces size
MEDIUM
