



## 1380

---

Throwaway disposable-fabric stretcher sheet - Green 230x90 cm



---

### FEATURES

Throwaway disposable-fabric stretcher sheet. Available in green (25 gr/mq - cm 90x230) or white (28 gr/mq - cm 100x220). 10 sheets per package.

---

### MAIN FEATURES

Width (cm)	90
Length (cm)	230
Material	TNT

---

### SEE ALSO



**1382**

Throwaway disposable-fabric  
stretcher sheet - White 220x100 cm

---