



1354

Red knapsack "115" with side pockets and compartments



FEATURES

Knapsack made of red anti-tear nylon, double seams, reflective strips, reinforced shoulder straps and strong 10 mm zips. Supplied with front opening wide space inside of which are 4 rack pockets to rationally store various devices. Possibility to fix 2 litres oxygen bottle. Two wide lateral sections contain four bags in different colours allowing fast identification of various kits. PVC Water proof bottom.
Dimensions (LxHxD): 44x56x26 cm

MAIN FEATURES

Width (cm)	44
Height (cm)	56
Depth (cm.)	26
Weight (kg)	4
Material	Nylon

IMAGES



SEE ALSO



1362
"Ambumed" rescue knapsack



1366
Red knapsack "115" with side
pockets and compartments
COMPLETE



Rescuemed 1343
Rescue knapsack

ACCESSORIES



1354-001
Carrying bag for knapsack "115"
