



1354

Red knapsack "115" with side pockets and compartments



FEATURES

The backpack contains all the items necessary for medication, resuscitation, and emergency therapy. Suitable for mountain communities, firefighters, forest bodies, and demanding excursions.

This medical material transport backpack has three separate bodies connected by hinges, equipped with a vertical lifting handle on the central body and carrying handles on the external bodies; inside its compartments, there are four envelopes to contain specific materials. Adjustable shoulder straps and abdominal belt in Cordura 1000; padding and containment pocket on the operator's back.

Dimensions (LxHxD): 36,6 x 46,4 x 17,6 cm

MAIN FEATURES

Width (cm)	36,6
Height (cm)	46,4
Depth (cm.)	17,6

IMAGES



SEE ALSO



Rescuemed 1343
Rescue knapsack



1362
"Ambumed" rescue knapsack



1366
Red knapsack "115" with side
pockets and compartments
COMPLETE
