



11064

Manual suction unit



FEATURES

HAND SUCTION UNIT

Designed to provide quick and effective pharyngeal and tracheal suction for emergency use.

The adjustable stroke knob provides different suction pressure.

Supplied with 250 ml bottle and 2 non sterile adult cannulas 17 cm (Ø 6 mm) and 19 cm (Ø 10 mm).

Max vacuum: -400 mmHg

Peak free airflow: >20 l/min

External diameter of catheter connection: 10 and 6 mm

Weight: 260 g

Size: 230x60x190 mm

MAIN FEATURES

Length (cm)	23
Height (cm)	19
Width (cm)	6
Weight (kg)	0,26

ACCESSORY OF



1366
Red knapsack "115" with side
pockets and compartments
COMPLETE
