



910/GIN PROOF Winner

Self loading trolley for incubators predisposed for certified fixing system



FEATURES

Designed to load and transport incubators on rescue vehicles. This trolley is part of "Winner" range, keeping the all same features. Made of anodized aluminium, it has same features as Art. 910/GIN with an additional predisposition for our certified fixing system Art. 955 Proof.

MAIN FEATURES

Pediatric	
Height (cm)	87,5
Max loading height (cm)	65
Min loading height (cm)	62
Folded height (cm)	26
Width (cm)	56
Length (cm)	194
Length with folded trolley (cm)	155
Weight (kg)	31
Material	Anodized aluminium
Loading capacity (kg)	170

ACCESSORIES



Rear fixing system

Rear fixing system for roll-in H=148 mm



Front blocking device

Front blocking device for automatic loading stretcher



12048

Certified interface for incubator fixation for self loading trolleys



12050

Extended trolley incubator interface
