



## 12054

---

Oximeb Maxi - Oxygen cylinder holder



---

### FEATURES

Oxygen cylinders holder made of strong, antitear nylon and inside padding to prevent damages from shocks. Supplied with a carrying handle and a side pocket to contain tube and mask. Suitable for 3 liters oxygen cylinders and, thanks to a handy fastening, it can be mounted on the side protection of the self loading stretchers and on the back rest of the sedan chairs.

---

### MAIN FEATURES

Length (cm)	43
Width (cm)	16
Depth (cm.)	12
Weight (kg)	0,3
Material	Nylon

---

### SEE ALSO



938

Oximeb oxygen cylinder holder

---

---

## ACCESSORY OF

---



**Mercury Lite 7082/4RG Proof**  
Certified aluminium variable height  
stretcher high version with 4  
rotating wheels



**Mercury Cinque 7075/4RG PROOF**  
Certified 5 levels self-loading  
stretcher with 4 swivel wheels



**Mercury Lite Cinque 7095/4RG  
Proof**  
Certified aluminium 5 levels self-  
loading stretcher with 4 swivel  
wheels



**Mercury Lite Cinque 7097/4RG  
Proof**  
Certified aluminium 5 levels self-  
loading stretcher with 4 swivel  
wheels - high version



**Mercury EL 7092 PROOF**  
Certified self-loading stretcher



**7265 Proof**  
Aluminium self loading stretcher  
"Frog lite Advance" with  
extractable handles and  
Trendelenburg