



12010

Phleboclysis stand for stretchers Mercury and Gecko



FEATURES

Telescopic phleboclysis stand for stretcher. Thanks to its fixing system it is possible to fix it in every point on the stretcher's frame. **Dimensions: minimum height: cm 65 maximum height: cm 110**

MAIN FEATURES

Min height (cm)	47,5
Max height (cm)	130,5
Weight (kg)	0,5
Material	Aluminium alloy

IMAGES



SEE ALSO



936/A

Phlebotomy stand for stretchers
Art. 1510 Proof and Art. 1515 Proof



936

Phlebotomy stand for automatic
loading stretcher Winner and Frog

ACCESSORY OF



Gecko 7122/AZ

"GECKO" self loading stretcher
entry-level light blue



Gecko 7122/V

"GECKO" self loading stretcher
entry-level green
