



## 12010

---

Phleboclysis stand for stretchers Mercury and Gecko



---

### FEATURES

Telescopic phleboclysis stand for stretcher. Thanks to its fixing system it is possible to fix it in every point on the stretcher's frame. **Dimensions: minimum height: cm 65 maximum height: cm 110**

---

### MAIN FEATURES

Min height (cm)	47,5
Max height (cm)	130,5
Weight (kg)	0,5
Material	Aluminium alloy

---

### IMAGES



---

## SEE ALSO

---



**936/A**

Phlebotomy stand for stretchers  
Art. 1510 Proof and Art. 1515 Proof



**936**

Phlebotomy stand for automatic  
loading stretcher Winner and Frog

---

---

## ACCESSORY OF

---



**Gecko 7122/AZ**

"GECKO" self loading stretcher  
entry-level light blue



**Gecko 7122/V**

"GECKO" self loading stretcher  
entry-level green

---