




## 9016 Galaxy

"GALAXY" Red spinal board




### FEATURES

Spinal board made of plastic material at high strength. Designed paying particular attention to the difficulties met during every day's interventions; they have 4 holes for the quick and total fixing of the head immobilizer and two cavities where the board lays on the floor, when the base is blocked in the traditional way, that allow to avoid damages to the rip-off straps during the usage of the spinal board or the accommodation in the ambulance. It has 12 handles for the transport and is supplied with 3 belts with rapid unhooking buckle (Art. 605/MEB). X-ray translucent.

Device certified according to the European harmonized safety standards UNI EN 1865.

### MAIN FEATURES

Size	Adults
Usage	Immobilization
Length (cm)	183
Width (cm)	41
Thickness (cm)	5
Weight (kg)	5,5
Loading capacity (kg)	150
Trauma	
Material	ABS

### SEE ALSO



**Galaxy 9012**  
"GALAXY" Orange spinal board



**Galaxy 9010**  
"GALAXY" Yellow spinal board

---

## ACCESSORIES

---



**Spider 690**

Belts immobilization system Spider  
for spinal board



**Jelly 629**

Orange standard head immobilizer  
for spinal board



**12144/N**

Carrying case for spinal board

---

---

## STANDARD ACCESSORIES

---



**605/MEB**

50 mm orange belt in one piece  
with plastic buckle and ME.BER.  
ribbon

---