



1362

“Ambumed” rescue knapsack



FEATURES

New modular knapsack, practical, compact, simple and suitable for rescue volunteers and professionals. Externally it is made of blue tear-resistant double sewing nylon, sturdy 10 mm zips, reflective high visibility stripes on all sides and a closed-cell PVC bottom, with plastic supports. On the inner side there is a protective hard shell structure made of thermoplastic material, which consists of a central wide space, characterized by a double-zip opening, and 4 external side pockets. The central space is supplied with Velcro® stripes for holding 4 coloured pockets, each one supplied with a large transparent window reinforced with double sewing, and 2 professional bags to store and organize different tools and materials. This knapsack is completed by a frontal pocket with inner straps to secure different objects, 3 anti-slip rubber handles (1 on the upper side and 2 on lateral sides) to carry by hand and, on the back side, 2 shoulder bands with reflecting stripes.

Ambumed is complete with 4 pockets (Art. 1362-002) and 2 professional bags (art. 1362-001).

Dimensions: 53 x 53 x 25 cm (LxHxD).

MAIN FEATURES

Width (cm)	53
Height (cm)	53
Depth (cm.)	25
Weight (kg)	4,1
Material	Nylon

IMAGES



SEE ALSO



1325/RB
Compact professional rescue
backpack Antartide



Rescuemed 1343
Rescue knapsack



1354
Red knapsack "115" with side
pockets and compartments

ACCESSORIES



1362-001/R
Red squared inside bag for
knapsacks



1362-001/G
Yellow squared inside bag for
knapsacks



1362-001/V
Green squared inside bag for
knapsacks



1362-002/A
Orange rectangular inside bag for
knapsacks



1362-002/B
Blue rectangular inside bag for
knapsacks



1362-002/R
Red rectangular inside bag for
knapsacks



1373
Green multifunctional medical roll
