



1365

Rescue bag Nevada



FEATURES

Bag made of cordura royal color. It is equipped with high visibility reflective stripes, sturdy double zips, a comfortable shoulder strap with metal hooks and handles to carry by hand. The bag is made up of a central space divided into 3 wide compartments with rigid adjustable dividers. Externally it is characterized by a frontal compartment with a pocket and fixed dividers, 3 side pockets with zip and two open compartments for the placement of bottles of various kinds. It is completed by a customizable name plate on the lateral side. The bottom is water-proof made of PVC with plastic support feet. The linings can be easily removed to allow internal cleaning and disinfection. **Dimensions: 50 x 25 x 28 cm (LxHxD).**

MAIN FEATURES

Length (cm)	50
Height (cm)	25
Depth (cm.)	28
Weight (kg)	1,5

SEE ALSO



1345
Empty emergency bag



1342
Emergency empty "Baby Bag"