



19152-130


"MTX-AR-130" Reduced height stretcher support with lateral movement and 130 mm ramp, without underneath space



FEATURES

"MTX-AR-130" Reduced height stretcher support with lateral movement and 130 mm ramp, without underneath space.

MAIN FEATURES

Side sliding	
Width (cm)	60,5
Length (cm)	208
Weight (kg)	42
Loading capacity (kg)	400
Plate height	11,3
Total height	14,3

ACCESSORIES



19123

Ribaltina 130 mm con attacchi per supporti e rullo



19115

Ribaltina da 180 mm con attacchi per supporti e rullo