



1468

Bottle for oxygen – 2 L



FEATURES

Bottles for oxygen therapy in conformity with F.U. norms, produced with stainless steel, refillable and complete of test certificate.
2 Lt. Bottle = 150 bar – 4,1 kg. – 300 Lt. gaseous – H 400 mm. – diam. 120 mm.

MAIN FEATURES

Height (cm)	40
Diameter	12

SEE ALSO



1470
Bottle for oxygen – 5 L



1475
Bottle for oxygen – 7 L



1465
Bottle for oxygen – 0,5 L
