



## 1475

---

Bottle for oxygen – 7 L



---

### FEATURES

Bottles for oxygen therapy in conformity with F.U. norms, produced with stainless steel, refillable and complete of test certificate.  
7 L Bottle = 200 bar – 10,8 kg. – 1400 Lt. gaseous – H 700 mm. – diam 140 mm.

---

### MAIN FEATURES

Height (cm)	70
Diameter	14

---

### SEE ALSO



**1468**  
Bottle for oxygen – 2 L



**1465**  
Bottle for oxygen – 0,5 L



**1478**  
Bottle for oxygen – 10 L

---