



## 1478

---

Bottle for oxygen – 10 L



---

### FEATURES

Bottles for oxygen therapy in conformity with F.U. norms, produced with stainless steel, refillable and complete of test certificate.  
10 L Bottle = 200 bar - 13,25 kg. – H 955 mm – Ø 137 mm - 2000 L gaseous

---

### MAIN FEATURES

Height (cm)	95,5
Diameter	13,7

---

### IMAGES



---

SEE ALSO

---



**1468**  
Bottle for oxygen – 2 L



**1465**  
Bottle for oxygen – 0,5 L



**1475**  
Bottle for oxygen – 7 L

---