



7015

Folding sky mattress black for bariatric MBP platform




FEATURES

Mattress for bariatric MBP platform made of expanded polyurethane and fireproof, waterproof and anti-tear sky cover. It is foldable into 2 parts for an easy storage and it is supplied with 6 Velcro[®] strips for anchoring it to the platform frame.

Dimensions: Length x Width x Thickness: 173 x 96 x 5,5 (max thickness 10) cm

*All Pictures shown are for **illustration purpose** only.

MAIN FEATURES

Bariatric	
Length (cm)	173
Width (cm)	96
Thickness (cm)	5,5 10

IMAGES



SEE ALSO



12080-001

Monting kit for bariatric platform
on self loading stretcher

ACCESSORY OF



12080

"MPB" wide bariatric platform for
self-loading stretchers Frog
Mercury and Gecko complete of
belts with reel and fixing for
stretcher



**MBS-01 MEBER BARIATRIC
STRETCHER 7270**

Self loading stretcher for bariatric
transport