



12080

"MPB" wide bariatric platform for self-loading stretchers Frog Mercury and Gecko complete of belts with reel and fixing for stretcher




FEATURES

The "MBP" platform is a system for Frog and Mercury self-loading stretchers that allows the safe transport of bariatric patients. The stainless steel frame along with the aluminium platform ensure lightness and sturdiness at the same time. The device is completed by belts with reel system on both sides and 4 fixing brackets for anchoring it to the stretcher frame.

Dimensions: (L x L x H) 1375 x 935 x 130 mm

Note: it is recommended to indicate the stretcher model on which the MPB platform will be installed, so to adjust the fixing clamps.

MAIN FEATURES

Bariatric	
Length (cm)	137,5
Width (cm)	93,5
Height (cm)	13
Weight (kg)	18
Loading capacity (kg)	350

IMAGES



SEE ALSO



**MBS-01 MEBER BARIATRIC
STRETCHER 7270**
Self loading stretcher for bariatric
transport



12080-001
Monting kit for bariatric platform
on self loading stretcher

ACCESSORIES



12080-001
Monting kit for bariatric platform
on self loading stretcher



7015
Folding sky mattress black for
bariatric MBP platform

ACCESSORY OF



MBS-01 MEBER BARIATRIC STRETCHER 7270
Self loading stretcher for bariatric transport



Mercury Cinque 7075/4RG PROOF
Certified 5 levels self-loading stretcher with 4 swivel wheels



Mercury Cinque 7077/4RG PROOF
Certified 5 levels self-loading stretcher with 4 swivel wheels - high version



Mercury Lite Cinque 7095/4RG Proof
Certified aluminium 5 levels self-loading stretcher with 4 swivel wheels



Mercury Lite Cinque 7097/4RG Proof
Certified aluminium 5 levels self-loading stretcher with 4 swivel wheels - high version



Mercury EL 7092 PROOF
Certified self-loading stretcher



Mercury EL 7094 PROOF
Certified self-loading stretcher - high version



7265 Proof
Aluminium self loading stretcher "Frog lite Advance" with extractable handles and Trendelenburg



Mercury 7070/4RG PROOF
Stainless steel variable height stretcher certified with 4 rotating wheels



Mercury 7072/4RG PROOF
Certified stainless steel variable height stretcher high version with 4 rotating wheels



Mercury Lite 7080/4RG Proof
Ambulance cot in aluminium with variable height and four swivel castors



Mercury Lite 7082/4RG Proof
Certified aluminium variable height stretcher high version with 4 rotating wheels



Frog 7210
"FROG" yellow self loading stretcher with trendelenburg



Frog 7210 PROOF
"FROG" yellow self loading stretcher with trendelenburg certified



Frog 7240 PROOF
"FROG PLUS" self loading stretcher with trendelenburg certified



Frog Lite 7260 PROOF
"FROG LITE" self loading aluminium stretcher with trendelenburg certified