



12080

"MPB" wide bariatric platform for self-loading stretchers Frog Mercury and Gecko complete of belts with reel and fixing for stretcher



FEATURES

The "MBP" platform is a system for Frog and Mercury self-loading stretchers that allows the safe transport of bariatric patients. The stainless steal frame along with the alluminium platform ensure lightness and sturdyness at the same time. The device is completed by belts with reel system on both sides and 4 fixing brackets for anchoring it to the stretcher frame. Dimensions: (L x L x H) 1375 x 935 x 130 mm

Note: it is recommended to indicate the stretcher model on which the MPB platform will be installed, so to adjust the fixing clamps.

MAIN FEATURES

Bariatric	₩
Length (cm)	137,5
Width (cm)	93,5
Height (cm)	13
Weight (kg)	18
Loading capacity (kg)	350

IMAGES



SEE ALSO



MBS-01 MEBER BARIATRIC STRETCHER 7270 Self loading stretcher for bariatric transport



12080-001 Monting kit for bariatric platform on self loading stretcher

ACCESSORIES



12080-001Monting kit for bariatric platform on self loading stretcher



7015Folding sky mattress black for bariatric MBP platform

ACCESSORY OF



MBS-01 MEBER BARIATRIC STRETCHER 7270

Self loading stretcher for bariatric transport



Mercury Cinque 7075/4RG PROOF Certified 5 levels self-loading stretcher with 4 swivel wheels



Mercury Cinque 7077/4RG PROOF

Certified 5 levels self-loading stretcher with 4 swivel wheels high version



Mercury Lite Cinque 7095/4RG Proof Certified aluminium 5 levels self-

loading stretcher with 4 swivel wheels



Mercury Lite Cinque 7097/4RG Proof

Certified aluminium 5 levels self-loading stretcher with 4 swivel wheels - high version



Mercury EL 7092 PROOF Certified self-loading stretcher



Mercury EL 7094 PROOF Certified self-loading stretcher high version

7265 Proof

Aluminium self loading stretcher "Frog lite Advance" with exctractable handles and Trendelenburg



Mercury 7070/4RG PROOF

Stainless steel variable height stretcher certified with 4 rotating wheels



Mercury 7072/4RG PROOF

Certified stainless steel variable height stretcher high version with 4 rotating wheels



Mercury Lite 7080/4RG Proof Ambulance cot in aluminium with variable height and four swivel castors



Mercury Lite 7082/4RG Proof Certified aluminium variable

height stretcher high version with 4 rotating wheels



Frog 7210

"FROG" yellow self loading stretcher with trendelenburg



Frog 7210 PROOF

"FROG" yellow self loading stretcher with trendelenburg certified



Frog 7240 PROOF
"FROG PLUS" self loading stretcher with trendelenburg certified



Frog Lite 7260 PROOF

"FROG LITE" self loading aluminium stretcher with trendelenburg certified