



## 11022/FA Cormos Plus

---

White complete medical suction unit with autoclavable canister



---

### FEATURES

The CORMOS PLUS suction units are easy, compact and battery powered medical devices. They differ from CORMOS series for their higher power erogated by a suction flow 30 liter/minute engine. The suction rate can be controlled thanks to a 800 mbar (80 kPa) manometer. It is supplied with: autoclavable canister, carrying MeBer green bag, 12V and 220V battery charger, suction probe with hose and wall rechargeable bracket with 10G capacity.

Item in accordance with ISO 10079-1; IEC 600601-1 3rd ed., approval e1 according to Directive 72/245/EEC, UNI EN 1789, IEC 60601-1-11. Classification according MDD 93/42/EEC: IIb.

Sizes 35x12x24 cm / Weight 4,6 kg / Autonomy on continuous functioning about > 45 minutes

---

### MAIN FEATURES

Width (cm)	35
Height (cm)	24
Depth (cm.)	12
Weight (kg)	4,6

---

## SEE ALSO

---



**Cormos Plus 11022/LN**  
White complete medical suction unit with canister for disposable unit



**Cormos 11010/LN**  
Portable aspirator for secretions with disposable liner

---

## ACCESSORIES

---



**11025**  
Serres disposable liners 1000 ml (x 12 pz./pcs.)



**11024**  
Serres disposable liners 1000 ml (x 6 pz./pcs.)



**11026**  
Serres disposable liners 1000 ml (x 36 pz./pcs.)



**11012**  
Autoclavable canister



**11028**  
Autoclavable canister



**11029**  
Wall charging/retention bracket for CORMOS PLUS suction units



**11032**  
Custodia verde per aspiratore CORMOS PLUS



**11034**  
Battery charger for Cormos 2 poles



**11023**  
10g certified wall bracket with integrated charging system for CORMOS and CORMOS PLUS

---

---

## STANDARD ACCESSORIES

---



**11014**

Suction rigid probe and hose for  
"CORMOS" and "CORMOS PLUS"  
aspirators



**11016**

Green carrying bag for suction unit

---