



11022/FA Cormos Plus

White complete medical suction unit with autoclavable canister



FEATURES

The CORMOS PLUS suction units are easy, compact and battery powered medical devices. They differ from CORMOS series for their higher power erogated by a suction flow 30 liter/minute engine. The suction rate can be controlled thanks to a 800 mbar (80 kPa) manometer. It is supplied with: autoclavable canister, carrying MeBer green bag, suction probe with hose and wall rechargeable bracket with 10G capacity.

Item in accordance with ISO 10079-1; IEC 600601-1 3rd ed., approval eI according to Directive 72/245/EEC, UNI EN 1789, IEC 60601-1-11. Classification according MDD 93/42/EEC: IIb.

Sizes 35x12x24 cm / Weight 4,6 kg / Autonomy on continuous functioning about > 45 minutes

MAIN FEATURES

Width (cm)	35
Height (cm)	24
Depth (cm.)	12
Weight (kg)	4,6

SEE ALSO



Cormos Plus 11022/LN
White complete medical suction
unit with canister for disposable
unit



Cormos 11010/LN
Portable aspirator for secretions
with disposable liner

ACCESSORIES



11025
Serres disposable liners 1000 ml (x
12 pz./pcs.)



11024
Serres disposable liners 1000 ml (x 6
pz./pcs.)



11026
Serres disposable liners 1000 ml (x
36 pz./pcs.)



11012
Autoclavable canister



11028
Autoclavable canister



11029
Wall charging/retention bracket for
CORMOS PLUS suction units



11032
Custodia verde per aspiratore
CORMOS PLUS



11034
Battery charger for CORMOS 2 poles



11023
10g certified wall bracket with
integrated charging system for
CORMOS and CORMOS PLUS

STANDARD ACCESSORIES



11014

Suction rigid probe and hose for
"CORMOS" and "CORMOS PLUS"
aspirators



11016

Green carrying bag for suction unit
