



## 11022/LN Cormos Plus

---

White complete medical suction unit with canister for disposable unit



---

### FEATURES

The CORMOS PLUS suction units are easy, compact and battery powered medical devices. They differ from CORMOS series for their higher power erogated by a suction flow 30 liter/minute engine. The suction rate can be controlled thanks to a 800 mbar (80 kPa) manometer. It is supplied with: 1000ml autoclavable canister for disposable liner (Serres), MeBer carrying bag, 12V and 220V battery charger, suction probe with hose and wall rechargeable bracket with 10G capacity.

Item in accordance with ISO 10079-1; IEC 600601-1 3rd ed., approval e1 according to Directive 72/245/EEC, UNI EN 1789, IEC 60601-1-11. Classification according MDD 93/42/EEC: IIb.

Sizes 35x12x24 cm / Weight 4,6 kg / Autonomy on continuous functioning about > 45 minutes

---

### MAIN FEATURES

Width (cm)	35
Height (cm)	24
Depth (cm.)	12
Weight (kg)	4,6

---

## SEE ALSO

---



**Cormos Plus 11022/FA**  
White complete medical suction  
unit with autoclavable canister



**Cormos 11010/LN**  
Portable aspirator for secretions  
with disposable liner

---

## ACCESSORIES

---



**11025**  
Serres disposable liners 1000 ml (x  
12 pz./pcs.)



**11024**  
Serres disposable liners 1000 ml (x 6  
pz./pcs.)



**11026**  
Serres disposable liners 1000 ml (x  
36 pz./pcs.)



**11028**  
Autoclavable canister



**11029**  
Wall charging/retention bracket for  
CORMOS PLUS suction units



**11032**  
Custodia verde per aspiratore  
CORMOS PLUS



**11034**  
Battery charger for Cormos 2 poles



**11023**  
10g certified wall bracket with  
integrated charging system for  
CORMOS and CORMOS PLUS

---

## STANDARD ACCESSORIES

---



**11014**  
Suction rigid probe and hose for  
"CORMOS" and "CORMOS PLUS"  
aspirators



**11016**  
Green carrying bag for suction unit