



11022/LN Cormos Plus

White complete medical suction unit with canister for disposable unit



FEATURES

The CORMOS PLUS suction units are easy, compact and battery powered medical devices. They differ from CORMOS series for their higher power erogated by a suction flow 30 liter/minute engine. The suction rate can be controlled thanks to a 800 mbar (80 kPa) manometer. It is supplied with: 1000ml autoclavable canister for disposable liner (Serres), MeBer carrying bag, suction probe with hose and wall recargeable bracket with 10G capacity.

Item in accordance with ISO 10079-1; IEC 600601-1 3rd ed., approval e1 according to Directive 72/245/EEC, UNI EN 1789, IEC 60601-1-11. Classification according MDD 93/42/EEC: IIb.

Sizes 35x12x24 cm / Weight 4,6 kg / Autonomy on continuous functioning about > 45 minutes

MAIN FEATURES

| Width (cm) | 35 |
|-------------|-----|
| Height (cm) | 24 |
| Depth (cm.) | 12 |
| Weight (kg) | 4,6 |

SEE ALSO



Cormos 11010/LN Portable aspirator for secretions with disposable liner



Cormos Plus 11022/FA White complete medical suction unit with autoclavable canister

ACCESSORIES



11162 Disposable 1-litre canister



11024 Serres disposable liners 1000 ml (x 6 pz/pcs.)



11025 Serres disposable liners 1000 ml (x 12 pz./pcs.)



11026 Serres disposable liners 1000 ml (x 36 pz/pcs.)



11012 Autoclavable canister



11013 Antibacterical filter



11034 Battery charger for Cormos 2 poles



11022-003 2 poles cable for aspirator CORMOS

STANDARD ACCESSORIES



11014
Suction rigid probe and hose for "CORMOS" and "CORMOS PLUS" aspirators



11032 Custodia verde per apiratore CORMOS PLUS



10g certified wall bracket with integrated charging system for CORMOS and CORMOS PLUS



11028 Autoclavable canister