




962

EN 1789 certified fixing system for stretcher/chair "EVOX"



FEATURES

Fixing system for self loading stretcher/chair EVOX art. 7112 Proof is entirely made of stainless steel to ensure corrosion resistance and high durability. Due to its sturdy structure is able to withstand shocks of up to 10G in all directions, as regulated by the UNI EN 1789 norm.



Device certified according to the European harmonized safety standards UNI EN 1789.

MAIN FEATURES

Length (cm)	161
Width (cm)	43,5
Height (cm)	13
Weight (kg)	6

ACCESSORY OF



Evox Proof 7112
"EVOX" certified chair cot with
mattress



Evox Proof 7112I
"EVOX" certified self-loading
stretcher convertible into chair
complete with mattress