



## 975

---

“Ergolift fixing system” - Fixing system for chair Extra Ergolift



---

### FEATURES

Fixing device designed to ensure the fixation of Extra Ergolift chair inside the ambulance patient compartment. It is made of stainless steel and it is corrosion resistant and easy to clean. The innovative monobloc configuration allows the operator to lock and unlock the chair fastly and easily.

The fixing system for Extra Ergolift chair can be also used as fixing system for the traditional Extra chair: 670/BR, 676 and 678.

Device certified according to the European harmonized standards UNI EN 1789.



---

### MAIN FEATURES

Width (cm)	58
Height (cm)	111
Depth (cm.)	23
Weight (kg)	8

---

## IMAGES

---



---

## SEE ALSO

---



**960**  
"Extra fixing system" - fixing system  
for chair Extra

---

## ACCESSORY OF

---



**Extra Ergolift 671**  
Powered stair chair



**Extra 670/BR**  
Track assisted stair chair with arms



**Extra Lite Standard 676**  
Foldable carrying chair with 4  
wheels, footrest and arms



**678 Extra Lite Basic**  
Extra Lite Basic, foldable carrying  
chair

---