



## 566

---

Black bag with 3 belts in 2 pieces made of black ribbon with orange MeBer logo for Extra chairs



---

### FEATURES

Black bag with 3 belts 607/N-MEB in 2 pieces made of black ribbon with orange MeBer logo for Extra chairs.

---

### ACCESSORIES



#### 607/N-MEB

Black belt in two pieces with plastic buckle and Me.Ber. band height 50mm

---

---

ACCESSORY OF

---



**Extra 670/BR**  
Track assisted stair chair with arms



**679**  
Narrow foldable sedan chair



**Extra Lite Standard 676**  
Foldable carrying chair with 4 wheels, footrest and arms



**678 Extra Lite Basic**  
Extra Lite Basic, foldable carrying chair



**Extra Ergolift 671**  
Powered stair chair

---