



1366

Red knapsack "115" with side pockets and compartments COMPLETE



FEATURES

The backpack contains all the items necessary for medication, resuscitation, and emergency therapy. Suitable for mountain communities, firefighters, forest bodies, and demanding excursions. This medical material transport backpack has three separate bodies connected by hinges, equipped with a vertical lifting handle on the central body and carrying handles on the external bodies; inside its compartments, there are four envelopes to contain specific materials. Adjustable shoulder straps and abdominal belt in Cordura 1000; padding and containment pocket on the operator's back.

Dimensions (LxHxD): 45,5 x 51,5 x 27 cm

In this version (1366), the knapsack is supplied with:

- 1 Adult silicone resuscitation bag
- 1 Silicone mask size 0
- 1 Silicone mask size 5
- 1 Silicone mask size 4
- 1 Silicone mask size 3
- 1 Silicone mask size 2
- 1 Silicone mask size 1
- 1 Silicone reservoir bag for adults/children
- 1 Guedel airway size 4, red cannula, 100 mm
- 1 Guedel airway size 3, yellow cannula, 90 mm
- 1 Guedel airway size 2, green cannula, 80 mm
- 1 Guedel airway size 1, white cannula, 70 mm
- 1 Guedel airway size 0, black cannula, 60 mm
- 1 Guedel airway size 00, blue cannula, 50 mm
- 1 Oxygen cylinder, 2 L, with pressure reducing valve
- 1 Adult oxygen mask with reservoir
- 1 POCKET MASK resuscitation mask
- 1 VACUUM PUMP pistol suction aspirator
- 2 500 ml saline solution bottles / CE
- 1 Lister scissors, 19 cm, DIN 58279-A190
- 1 Adjustable adult cervical collar
- 1 Adjustable pediatric cervical collar
- 2 Burn dressings 40 x 60 cm, DIN 13152-BR
- 5 Pairs of sterile copolymer gloves, one size
- 2 "VISUAL" splash protection visors
- 2 Medical waste bags
- 1 125 ml disinfectant, 10% iodopovidone iodine, PMC
- 2 Packs of 12 sterile gauzes, 18 x 40 cm
- 10 Packs of 25 sterile gauzes, 10 x 10 cm
- 1 Case with 20 assorted plasters
- 1 Silk adhesive tape roll, 5 m x 2.5 cm
- 2 Hemmed bandages, 3.5 m x 10 cm
- 1 Sterile tweezers, 10 cm
- 1 Sterile drape, 60 x 80 cm, DIN 13152-A
- 2 Gold/silver thermal blankets, 160 x 210 cm
- 2 Instant ice packs
- 1 Digital thermometer
- 1 PERSONAL sphygmomanometer with stethoscope
- 2 EMOCONTROL anti-hemorrhagic bandages
- 3 Latex-free tourniquets

MAIN FEATURES

Width (cm)	45,5
Height (cm)	51,5
Depth (cm.)	27

IMAGES



SEE ALSO



1354
Red knapsack "115" with side
pockets and compartments
